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TEL06-6110-4100

# Medi·Ca CC

## Ready-made dry medium for coliform count

### Introduction

Medi-Ca is a microbiological culture device made up of a waterproof sheet, a ready-made dry medium on the sheet and a transparent cover over the medium. The Medi-Ca CC method is intended to indicate the level of coliform bacteria in food products as an alternative to the Violet Red Bile Agar method in U.S. Food and Drug Administration's Bacteriological Analytical Manual (FDA-BAM): Enumeration of *Escherichia coli* and the Coliform Bacteria. It is compact and easy to use, and reduces the total amount of waste produced during testing. The Medi-Ca CC method has been granted AOAC Research Institute *Performance Tested Methods*<sup>SM</sup> (021401).

### Procedure

#### 1. Sample preparation

- 1) Weigh each 50 g test portion into a blender jar.
- 2) Add 450 ml Butterfield's phosphate-buffered diluent (BPD) and blend for up to 2 minutes. For a recipe of BPD, see the FDA website (BAM Reagents R11).
- 3) Prepare all decimal dilutions with 90 ml BPD plus 10 ml previous dilution.
- 4) Shake 25 times in a 30 cm arc.

#### 2. Inoculation

- 1) Remove each Medi-Ca CC sheet from an aluminum bag under aseptic conditions.
- 2) Place the sheet on a flat surface and allow it to reach room temperature (15-25°C).
- 3) Lift the cover and place 1 ml sample suspension onto the center of the medium.
- 4) Drop the cover onto the sample.

**Do not use the sheet with a bent cover to ensure the sample spreads on the entire medium.**

- 5) Leave the sheet for 3 minutes or more on a horizontal surface.

**Do not move or tilt the sheet until solidification of the suspension is completed to avoid the sample being spilled out from the growth area.**

- 6) Hold both ends of the sheet and place it into an incubator.

**The sheets can be held after the solidification without any time restriction.**

#### 3. Incubation

Incubate the sheets at  $35 \pm 1^\circ\text{C}$  for  $24 \pm 1$  hours. It is possible to stack up to 25 sheets.

#### 4. Interpretation

Count all blue colonies regardless of size or intensity. The suitable colony counting range is 1-250.

#### Trouble shooting for counting colonies:

- a) When a number of colonies per sheet exceeds 250, for all dilutions, record the count as too numerous to count (TNTC). If an estimated count is required, count colonies within 1-3 squares (1 cm × 1 cm) printed on the cover and calculate an average. Multiplying the average number by 20 provides the estimated count since the circular growth area is approximately 20 cm<sup>2</sup>.
- b) When the entire growth area becomes blue, record the count as too numerous to count (TNTC).
- c) When a bubble disrupts a colony so that the colony outlines the bubble, count it as one colony.
- d) When a colony is spreading, count it as one colony.
- e) When two or more spreading colonies appear to originate from separate sources, count each source as one colony.

- f) When the sample is not clear (i.e. cloudy or dark), prepare a higher dilution.
- g) When the entire growth area become blue due to food components involving the chromogenic reaction, prepare a higher dilution.
- h) When colonies are extremely small or light, incubate continuously for a few hours to make the colony size larger or the colony color intensity darker.

#### 5. Colony isolation

Lift the cover and pick a single colony from the gel.

### Cautions

- 1) This product is not for clinical testing.
- 2) This product has been validated for raw meat, raw poultry, raw fish and fresh fruits in AOAC *Performance Tested Methods*<sup>SM</sup> program, but has not been tested with all possible food products, food processes, test protocols or strains.
- 3) Do not open the cover until just before inoculation.
- 4) Do not use the product after its expiry date.
- 5) Do not use the product that is damaged, deformed or discolored, or that contains any foreign materials.
- 6) Do not expose the product to direct sunlight.
- 7) Do not push the cover after inoculation to avoid sample suspension being spilled out from a growth area.
- 8) If sample suspension spilled out from a growth area, try again with a new sheet.
- 9) If medium or reagent gets into eyes or mouth, rinse immediately with plenty of water and consult a doctor.
- 10) Analysis needs to be performed under a laboratory analyst with microbiological training and supervision. All waste must be handled as a biohazard and disposed by autoclaving.

### Storage

Keep in the refrigerator (2-8°C).

### Storage after opening the bag

Put any unused sheet back into the bag, fold the end of the bag over twice, and seal with tape. The shelf life under refrigerated condition is 3 months after opening.

### Shelf life

Expiry date is specified on the product label after the word "EXP". The shelf life of the product is 18 months after manufacturing if it is stored properly without opening the bag.

### Disposal

Any and all media, supplements, and reagents must be sterilized by autoclaving after use, and then disposed as industrial waste according to local laws and regulations.

### Disclaimer of warranties

If any Medi-Ca CC sheet is proved to be defective by reason of improper workmanship and/or material, they will be replaced with new ones. It is the user's responsibility to determine that any test method and result meet its customers' or suppliers' requirements.

### Contact information

Business Development Department, Global Business Division,  
Packaging Operations, Dai Nippon Printing Co., Ltd.

Address: 1-1-1, Ichigaya-Kagacho, Shinjuku-ku, Tokyo 162-8001, Japan

URL: [https://www.dnp.co.jp/CGL/inquiry\\_eng/form.cgi](https://www.dnp.co.jp/CGL/inquiry_eng/form.cgi)

**Manufactured by**

**Dai Nippon Printing Co., Ltd.**

**1-1-1, Ichigaya Kagacho, Shinjuku-ku, Tokyo, 162-8001, Japan**